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1. Commercial name:

EasyClip
 EasyClip+
 Lotus
 Self-ligating single buccal tubes

2.  Manufacturer:

Aditek do Brasil Ltda.
 Address: Cesário Mota, 14
 ZIP: 14140-000 - Cravinhos – SP - Brazil
 Phone: +55 (16) 3951-9355 / 0800-0159520
 CNPJ 64.602.097/0001-95
 Technical Responsible: José Arantes Nogueira Neto - CREA: 601037885
 www.aditek.com.br
 aditek@aditek.com.br

3.  1304 Notified Body:

SIQ Ljubljana
 SLO-1000 Ljubljana, Trzaska cesta 2
 tel.: +386 1 4778 100
info@siq.si
www.siq.si

4. Intended Use:

The use of orthodontic Brackets is restricted to dental professionals with qualification in orthodontics. Use by persons outside this context may result in irreversible damage to the patient.

Aditek products have been developed for use in the oral cavity environment, directly on the teeth, with the objective of supporting orthodontic treatments in humans, promoting a better positioning of the teeth in the dental arches toward a correct position and occlusion.

The brackets are used in conjunction with the archwires, wires, springs, bands, elastomers and other auxiliary products for the orthodontic treatment.

The self-ligating brackets and tubes have a clip that holds the wire into the slot, eliminating the use of ligatures or ligature wires.

They are indicated for treating patients from 6 years old who have malocclusion problems such as prognathism, retrognathism, deep bite, open bite, crowding, anterior crossbite, posterior crossbite and diastema.

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NOTE: The choice of material to be used is made by the qualified professional, based on his knowledge acquired during his training in conjunction with the diagnosis and treatment plan for each patient.

5. Users:

The application of orthodontic Brackets is restricted to dental professionals with qualification in Orthodontics. Use by persons outside this context may result in irreversible damage to the patient.

6. Performance Characteristics:

Commonly orthodontic treatment without interruption takes from 1 to 3 years, varying in each case according to the facial pattern, type of malocclusion and oral conditions of the patient, which can be extended according to the complexity of the case.

7. Clinical Benefits:

The main benefits that these systems bring to the life of the patient is the aesthetic improvement with the alignment and leveling of the teeth, better hygiene, more efficient masticatory function facilitating swallowing, better digestion, prevention of occlusal problems, better speech and breathing.

8. Warning:

Product containing Chromium and Nickel in its composition, and may cause allergic reactions in patients susceptible to allergy to these substances. In case of reaction, remove the appliance and instruct the patient to seek medical care.

Aditek's products are designed and manufactured for single use only. The reuse or reprocessing may cause cross infection, decreased efficiency in bonding and loss of mechanical properties due to natural wear and forces applied in the removal of the product.

If the package has been breached, it is recommended to discard the product, minimizing the risk of transmission of microorganisms or possible product adulteration.

Store product in cool, dry and out of direct sunlight, keeping it in its own package until use.

The product is not recommended in the MRI environment, since the product will reflect the magnetic waves generated by the tomograph, preventing or masking the correct reading of the examination, so they should always be removed before MRI of the head and neck.

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The product is safe for use in pregnant women. It is known that due to hormonal changes during pregnancy the woman becomes more susceptible to having periodontal problems. During orthodontic treatment, the difficulty of cleaning can aggravate these problems. It is important that the patient is referred for professional cleaning at regular intervals.

The procedure for bonding and removing the brackets from the teeth must be performed according to the instructions of the manufacturer of the acquired orthodontic adhesive.

The composition of the product does not require specific procedures for transport and presents no risk to people.

The product is supplied in the non-sterile condition, however it has proven biocompatibility and no cleaning, disinfection or sterilization measures are necessary.

The professional should use good clinical practices for the disposal of medical products in accordance with current sanitary regulations.

Aditek do Brasil Ltda. is not responsible for the results obtained by using the product without necessary precautions or non-compliance with the warnings and instructions for use.

In the event of serious incidents with the product, Aditek do Brasil must be notified to the national competent authority in which the orthodontist is established.

9. Contraindications:

- Patients with poor oral hygiene,
- Patients with inability to cooperate with the treatment,
- Patients who have allergy to chromium and / or nickel,
- Disease or limitation that might hinder treatment,
- Existing bone or root resorption,
- Existing decalcification of dental enamel,
- Existing periodontal complications.

It is the primary responsibility of the orthodontist to identify any possible contraindication that may preclude the use of the product.

10. Residual Risks

All finished orthodontic treatments are subject to relapse. After treatment, the control must be determined by the professional to avoid or minimize relapse.

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11. Side Effects

The use of the medical device may result in discoloration or decalcification of teeth, bone and root resorption, periodontal complications, allergic reactions, difficulties in maintaining oral hygiene, oral and mucosal damage, difficulty in speaking or chewing, discomfort and pain.

12. Composition:

Stainless Steel AISI 17-4 – Bracket								
Limits of Chemical elements (%)								
C	Mn	P	S	Si	Cr	Ni	Cu	Fe
≤ 0,070	≤ 1,00	≤ 0,040	≤ 0,030	≤ 1,00	15,0 a 17,5	3,0 a 5,0	3,0 a 5,0	Balance
ASTM: A564								

Níquel Titânio NiTi #1 - Clip									
Limits of Chemical elements (%)									
Ni	C	Co	Cu	H	Fe	Nb	N + O	Cr	Ti
54,5 a 57,0	≤ 0,050	≤ 0,050	≤ 0,010	≤ 0,005	≤ 0,050	≤ 0,025	≤ 0,050	≤ 0,010	Balance
ASTM: F2063									

13. Shelf life:

10 years

14. Symbols and Meanings:

	Manufacturer		Authorized representative in the European Community
	Consult instructions for use		Complies with the European directive. XXXX represents the number of the certifying body
	Attention		Product contains chrome and nickel
	Do not use if the package is damaged or breached		Do not reuse product
	Product Reference Code		Lot number
	Expiration date		

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15. Patient information:

Product containing Chromium and Nickel in its composition, which may cause allergic reactions in people susceptible to these components;

Perform daily good hygiene to prevent plaque installation and subsequent periodontal disease;

It is important that the pregnant patient be referred for professional cleaning at regular intervals, due to the increased susceptibility of periodontal problems during pregnancy;

Chew hard food can cause the devices break or fall;

In case of damage, look for an orthodontist;

If swallowed or aspirated, look for medical care;

When participating in sports, always wear mouth guards, as impacts can damage appliances and harm in the oral cavity;

Mention all the undesirable side effects;

Mention about the risk of recurrence after treatment;

The patient must follow the professional's instructions regarding retainer to avoid recurrence;

The product is not recommended in the MRI environment.

16. Instructions for use:

The steps below must be performed only by dental professionals with qualification in Orthodontics.

Prepare the tooth surface according to the instructions of the manufacturer of the acquired orthodontic adhesive.

Bond the brackets on the teeth according to the instructions of the manufacturer of the acquired orthodontic adhesive.

Insert the instrument (explorer) in the circular opening of the clip located in between the lower tie-wings perpendicular to the vestibular face of the tooth.

Gently rotate the instrument to the left or to the right; in case of brackets with hooks, rotate the instrument in the opposite direction of the hook.

The clip may also be open using a fine tip explorer with a gentle lever action downward. The lower bracket clips open to gingival and the upper to occlusal.

Engage the selected archwire.

Place one tip of the tweezers over the upper tie-wings and the other tip below the clip.

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Press down the clip vertically. The first "click" unlocks the clip. The second, stronger click will lock the clip, closing the slot.

After treatment, remove the brackets and orthodontic adhesive from the teeth using an appropriate instrument as per the orthodontic adhesive manufacturer's instructions for use.

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