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www.siq.si

4. Intended Use:

The use of Orthodontic Bands is restricted to dental professionals with qualification in orthodontics. Use by persons outside this context may result in irreversible damage to the patient.

Aditek products have been developed for use in the oral cavity environment, directly on the teeth, with the objective of supporting orthodontic treatments in humans, promoting a better positioning of the teeth in the dental arches toward a correct position and occlusion.

The Bands are used in conjunction with the brackets, archwires, wires, springs, bands, elastomers and other auxiliary products for the orthodontic treatment.

They are indicated for treating patients from 6 years old who have malocclusion problems such as prognathism, retrognathism, deep bite, open bite, crowding, anterior crossbite, posterior crossbite and diastema.

NOTE: The choice of material to be used is made by the qualified professional, based on his knowledge acquired during his training in conjunction with the diagnosis and treatment plan for each patient.

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5. Users:

The application of Orthodontic Bands is restricted to dental professionals with qualification in Orthodontics. Use by persons outside this context may result in irreversible damage to the patient.

6. Performance Characteristics:

Commonly orthodontic treatment without interruption takes from 1 to 3 years, varying in each case according to the facial pattern, type of malocclusion and oral conditions of the patient, which can be extended according to the complexity of the case.

7. Clinical Benefits:

The main benefits that these systems bring to the life of the patient is the aesthetic improvement with the alignment and leveling of the teeth, better hygiene, more efficient masticatory function facilitating swallowing, better digestion, prevention of occlusal problems, better speech and breathing.

8. Warning:

Product containing Chromium and Nickel in its composition, and may cause allergic reactions in patients susceptible to allergy to these substances. In case of reaction, remove the appliance and instruct the patient to seek medical care.

Aditek's products are designed and manufactured for single use only. The reuse or reprocessing may cause cross infection, decreased efficiency in bonding and loss of mechanical properties due to natural wear and forces applied in the removal of the product.

If the package has been breached, it is recommended to discard the product, minimizing the risk of transmission of microorganisms or possible product adulteration.

Store product in cool, dry and out of direct sunlight, keeping it in its own package until use.

The product is not recommended in the MRI environment, since the product will reflect the magnetic waves generated by the tomograph, preventing or masking the correct reading of the examination, so they should always be removed before MRI of the head and neck.

The product is safe for use in pregnant women. It is known that due to hormonal changes during pregnancy the woman becomes more susceptible to having periodontal problems. During orthodontic treatment, the difficulty of cleaning can aggravate these problems. It is important that the patient is referred for professional cleaning at regular intervals.

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The procedure for bonding and removing the brackets from the teeth must be performed according to the instructions of the manufacturer of the acquired orthodontic adhesive.

The composition of the product does not require specific procedures for transport and presents no risk to people.

The product is supplied in the non-sterile condition, however it has proven biocompatibility and no cleaning, disinfection or sterilization measures are necessary.

The professional should use good clinical practices for the disposal of medical products in accordance with current sanitary regulations.

Aditek do Brasil Ltda. is not responsible for the results obtained by using the product without necessary precautions or non-compliance with the warnings and instructions for use.

In the event of serious incidents with the product, Aditek do Brasil must be notified to the national competent authority in which the orthodontist is established.

9. Contraindications:

- Patients with poor oral hygiene,
- Patients with inability to cooperate with the treatment,
- Patients who have allergy to chromium and / or nickel,
- Disease or limitation that might hinder treatment,
- Existing bone or root resorption,
- Existing decalcification of dental enamel,
- Existing periodontal complications.

It is the primary responsibility of the orthodontist to identify any possible contraindication that may preclude the use of the product.

10. Residual Risks

All finished orthodontic treatments are subject to relapse. After treatment, the control must be determined by the professional to avoid or minimize relapse.

11. Side Effects

The use of the medical device may result in discoloration or decalcification of teeth, bone and root resorption, periodontal complications, allergic reactions, difficulties in maintaining oral hygiene, oral and mucosal damage, difficulty in speaking or chewing, discomfort and pain.

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 10 years

14. Symbols and Meanings:

| | | | |
|---|--|---|---|
|  | Manufacturer |  | Authorized representative in the European Community |
|  | Consult instructions for use |  | Complies with the European directive. XXXX represents the number of the certifying body |
|  | Attention |  | Product contains chrome and nickel |
|  | Do not use if the package is damaged or breached |  | Do not reuse product |
|  | Product Reference Code |  | Lot number |
|  | Expiration date | | |

15. Patient information:

Product containing Chromium and Nickel in its composition, which may cause allergic reactions in people susceptible to these components;

Perform daily good hygiene to prevent plaque installation and subsequent periodontal disease;

It is important that the pregnant patient be referred for professional cleaning at regular intervals, due to the increased susceptibility of periodontal problems during pregnancy;

Chew hard food can cause the devices break or fall;

In case of damage, look for an orthodontist;

If swallowed or aspirated, look for medical care;

When participating in sports, always wear mouth guards, as impacts can damage appliances and harm in the oral cavity;

Mention all the undesirable side effects;

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Mention about the risk of recurrence after treatment;

The patient must follow the professional's instructions regarding retainer to avoid recurrence;

The product is not recommended in the MRI environment.

16. Instructions for use:

The steps below must be performed only by dental professionals with qualification in Orthodontics.

- Select the size of the band, choosing the one that best fits the tooth.
- Weld the accessories to be used.
- Prepare the tooth surface according to the instructions of the manufacturer of the acquired orthodontic adhesive or cement.
- Apply the resin or orthodontic cement inside the band, according to the manufacturer's instructions for use of the adhesive or cement purchased,
- Place the band on the teeth using specific instruments, according to the manufacturer's instructions for use of the adhesive or cement purchased.

After treatment, remove the band and orthodontic adhesive from the teeth using an appropriate instrument as per the orthodontic adhesive manufacturer's instructions for use.

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